Exemplar 2pCoin Analysis

Summary of Results from experiment

These results are heavily simplified JUST to show the values required for analysis of the experiment.

Coin Thickness	1	2	3	4	5	Ave
Reading /mm	1.84	1.86	1.81	1.80	1.85	1.83
E coin Diameter - 126mm - Coin diameter - 126/5 - 25.2mm						

5 coin Diameter = $126mm \rightarrow$ Coin diameter = 126/5 = 25.2mm



	Density / kgm ⁻³
1 coin	7320
From 5 coins	7540

Mass	Х/	Υ/	M _{coins} /g
Calculation	mm	mm	
Single coin	206	135	6.53
Stack of 5	208	143	34.4
coins			(1coin = 6.88g)

Analysis

(This is an example! Your values will be different BUT the methods should be the same)

State the (absolute) uncertainty in your value for the thickness of the coin, h.

Would either be the $\frac{1}{2}$ smallest division of the micrometer at either end (there is also an uncertainty in where the micrometer closes, so $2 \times \frac{1}{2}$ of 0.01mm i.e. 0.01 mm) OR $\frac{1}{2}$ range of repeats... In this case the <u>range</u> of readings of coin thickness is 1.86-1.80 = 0.06mm. The $\frac{1}{2}$ range = 0.03mm so $\frac{1}{2}$ range is greater so <u>Absolute uncertainty of coin thickness = 0.03mm</u>

Calculate the percentage experimental uncertainty in your value for h using: h = 1.83mm; absolute uncertainty = $0.03mm \rightarrow \%U = (0.03/1.83)*100 = 1.64\%$

In measuring the diameter of the 2p coins, explain why it is better to measure across all of the coins rather than measuring them individually

The 5 coin measurement is larger and therefore will have a lower % uncertainty. It will also reduce the effect of any random error associated with differences in the diameter of individual coins.

State the (absolute) uncertainty in your value for the measurement of 5d. The ruler's smallest division is 1mm ($\frac{1}{2}$ is 0.5mm) BUT counted at the start and end so doubled so uncertainty of diameter incorporating 2 measurements = 1mm

Calculate the percentage experimental uncertainty in your value for 5d: 5d = 126mm; absolute uncertainty = $1mm \rightarrow \%U = (1/126)*100 = 0.794\%$

Work out the total percentage experimental uncertainty for the volume value you calculated. (% uncertainty in volume = $(2 \times \%$ uncertainty in d) + % uncertainty in h) $\%U_h=1.64\%$; $\%U_d=1.59\%$; $\Rightarrow \underline{\%U_{vol}} = (2 \times 0.794) + 1.64 = 3.23\%$

State the (absolute) uncertainty in your value for distance "x" The ruler's smallest division is 1mm so the uncertainty of "x" incorporating 2 measurements = $\underline{1mm}$

Calculate the percentage uncertainty in your value for distance "x" for 1 coin and 10g x = 206mm; absolute uncertainty = 1mm \rightarrow %U =(1/206)*100 = 0.485%

for 5 coins and 50g

x = 208mm; absolute uncertain	ty = 1mm	→ %U =	(1/208)*100 =	0.481%
-------------------------------	----------	--------	--------	---------	--------

State the (absolute) uncertainty in your value for distance "y"

The ruler's smallest division is 1mm so the uncertainty of "y" incorporating 2 measurements = 1mm

Calculate the percentage uncertainty in your value for distance "y" for 1 coin and 10g

y = 135mm; absolute uncertainty = 1mm \rightarrow %U =(1/135)*100 = 0.741%

for 5 coin and 50g

y = 143*mm*; absolute uncertainty = 1*mm* → %*U* =(1/143)*100 = 0.699%

Work out the total percentage experimental uncertainty for the mass value you calculated

(if we assume %U_{Mass(mass)} =0%) then %U_{mass(coin)} = %U_x + %U_y

for 1 coin and 10g $%U_{mass(coin)} = %U_x + %U_y = 0.485 + 0.741 = 1.23 \%$

for 5 coin and 50g $%U_{mass(coin)} = \%U_x + \%U_y = 0.481 + 0.699 = 1.18\%$

Calculate the % difference between the two values of coin mass you obtained in the experiment.

% difference = <u>difference between the 2 values</u> × 100 average of the 2 values

mass of coin (from 1 coin and 10g) = 6.53g; mass of coin (from 5 coin and 50g) = 6.88g Average = (6.53+6.88)/2 = 6.71g; Difference = 6.88-6.53 = 0.35g% diff = (0.35/6.71)*100 = 5.22% (that means they were consistent with each other)

Work out the total percentage experimental uncertainty for the density value you calculated in the previous experiment $U_{density} = U_{mass} + U_{volume}$

For 1 coin and 10g (NB %U_{volume} is the same for both!) % $U_{mass} = 1.48 \%$ % $U_{volume} = 3.23\%$ % $U_{density} = \%U_{mass} + \%U_{volume} = 1.23 + 3.23 = 4.46 \%$

For 5 coins and 50g

 $U_{mass} = 0.955 \%$ $U_{volume} = 3.23\%$ $U_{density} = U_{mass} + U_{volume} = 1.18 + 3.23 = 4.41\%$

Work out the maximum and minimum values for your density value. (This is just going to be done for the 5coin and 50g value) $Density=7540 \text{ kgm}^{-3}$ Maximum= $7540 \times 1.044 =$ <u>7870 kgm}^{-3}</u> Minimum = $7540 \times 0.956 =$ <u>7210 kgm</u>⁻³

Which metal is it most likely to be?

The only metal from the table that has a density that overlaps with this result is Mild Steel

Metal	Density/kgm ⁻³
Copper	8930
Brass	8500
Mild steel	7860